



Icahn  
School of  
Medicine at  
**Mount  
Sinai**

*The Tisch  
Cancer Institute*



*Black Family  
Stem Cell Institute*

*Skin Biology  
and Diseases  
Resource-based  
Center*


*Department of Cell,  
Developmental and  
Regenerative Biology*




# EMBL epigenetics Practicum Introduction

---



# Melanoma is a malignancy of the pigment-producing cells melanocytes



# ARID2 is one of the most frequently mutated genes in melanoma



# ARID2 mutations spectrum in melanoma




## ARID2

RefSeq: [NM\\_152641](#)  
Ensembl: [ENST00000334344](#)  
CCDS: [CCDS31783](#)  
UniProt: [ARID2\\_HUMAN](#)


Somatic Mutation Frequency: 16.8% [i](#)

58 Missense 46 Truncating  
0 Inframe 1 Other

# The SWI/SNF complexes their and ATPase subunits



# SWI/SNF complex subunits are mutated in multiple cancers



ARID2 loss alters SWI/SNF occupancy and chromatin accessibility which alters transcription factor binding patterns, and results in aberrant transcriptional programs leading to increased metastases.




# Model

- Cell line:
  - SKmel147 – melanoma cell line
  - NRAS Q61R

- Conditions:
  - WT (CRISPR sgControl)
  - KO (CRISPR sgARID2)

- Growing conditions:
  - 37 °C, 5% CO<sub>2</sub>
  - DMEM (Dulbecco Media), %5 penicillin/streptomycin, 10% Fetal bovine serum

- Assay
  - ATAC-seq



\*Samples used for ATAC-seq analysis



# Assignment expectations

---

- Perform quality control, alignment and filtering of data generated from ATAC-seq
- Perform differential chromatin accessibility analysis and integrate with RNAseq differential expression results.
- Visualize differential chromatin accessibility results and its integration with RNAseq results.
- Perform motif analysis of differentially accessible regions.
- Perform functional analysis of genes associated with differentially accessible regions
- Annotate chromatin accessible regions with their functional genomic region and their nearest gene (Homer)

# Module goals

---




- Develop a basic understanding on best practices for quality assessment of chromatin profiling methods.
- Become familiar with the processing steps of ATAC-seq analysis, the aim of each step and the specific arguments/functions used.
- Obtain a basic understanding of the biology of chromatin remodelers, transcription factor binding, chromatin profiling and their dynamic correlation.
- Familiarization with basic visualization practices for chromatin profiling methods.

# Questions to drive the analysis

---

- Characterize the chromatin landscape at open regions?
- Which transcription factors are enriched at open regions? How many differentially accessible regions are between ARID2 WT and KO cells?
- What genes/pathways are enriched nearby the differential accessible regions?
- Which transcription factors are enriched at differential regions?
- What is the relationship between differentially open regions and differentially expressed genes?
- Hypothesize about the biological output (phenotype) produced by that changes that occur at the chromatin level.

# Workflow



# Pre-processing ATACseq Steps

---

- Adaptor removal: **TrimGalore**
- Alignment to GRCH38: **Bowtie2**
- Filtering/Sorting/Indexing sam files: **Samtools**
- Duplicate Removal: **Picard**
- Sorting/Indexing sam files: **Samtools**
- Generation of coverage tracks: **Deeptools**
- Calling significant peaks: **Macs2**
- Heatmap Visualization: **Deeptools**
- Coverage Visualization: **UCSC browser/IGV**
- Quantification of reads in peak: **Subread**
- Differential peak analysis: **AWK/other**
- Motif enrichment analysis: **Homer**
- Gene peak association (RNAseq results): **Bedtools**
- Pathway analysis: **Enrichr**

# Dataset integration (RNAseq/ChIPseq/ATACseq)

---



- Association of differential ATACseq regions with differential genes
- Visualization of ChIPseq enrichment at differential regions

# Files provided – Data to analyze

---



- Fastq files (ATACseq)
  - SKmel147 ARID2 WT ATAC-seq (2 files R1 and R2)
  - SKmel147 ARID2 KO ATAC-seq (2 files R1 and R2)
- Bigwig files (ChIPseq)
  - SKmel147 ARID2 WT ARID2
  - SKmel147 ARID2 WT FOSL2
  - SKmel147 ARID2 WT H3K4me3
  - SKmel147 ARID2 WT H3K27ac
- Table
  - RNAseq results for the differential gene expression analysis between SKmel147 ARID2 WT vs SKmel147 ARID2 KO



# Files provided – Containers and software

---



## #Containers

- bedtools\_v2.31.1
- bowtie2\_v2.4.5.sif
- deeptools\_v3.5.1.sif
- fastqc\_v0.11.8.sif
- homer.sif
- macs\_v2.2.9.1.sif
- picard\_v2.9.2.sif
- samtools\_v1.15.sif
- subread\_v2.0.1.sif
- trim-galore\_v0.6.9.sif

## #Other

- Singularity/3.11.0
- AWK
- IGV

# Files provided – Supporting files

---



- GRCH38 blacklisted regions
  - Regions defined by the Encode consortium with challenging mapping characteristics.
- GRCH38 -\+ 1 Kb promoters
  - Bed files of the promoter of genes for the GRCH38 genome version
- Bowtie2 Index
  - Files need for alignment of FASTQs with Bowtie2 corresponding to the GRCH38 genome version
- Code
  - “Draft” version of code to be completed with the resources provided for the course.

# Files – Paired end ATAC-seq

## ATAC-seq

- Melanoma cell line SKmel147
- Reads only from chromosome 5
- CRISPR ARID2KO vs ARID2 WT cell lines

File name	Cell line	Genotype
SKMel147_ARID2WT_ATAC_chr5_R1.sam	SKmel147 (Melanoma)	ARID2 WT
SKMel147_ARID2KO_ATAC_chr5_R1.sam	SKmel147 (Melanoma)	ARID2 KO

# Files – Single end ChIPseq

## ChIPseq-seq

- Melanoma cell line SKmel147
- Reads only from chromosome 5
- Parental cell lines

## Bigwigs

File name	Cell line	Genotype
SKmel147-H3K4me3_chr5.bw	SKmel147 (Melanoma)	ARID2 WT
SKmel147-H3K27ac_chr5.bw	SKmel147 (Melanoma)	ARID2 WT
SKmel147-FOSL2_chr5.bw	SKmel147 (Melanoma)	ARID2 WT
SKmel147-ARID2_chr5.bw	SKmel147 (Melanoma)	ARID2 WT

## ChIPseq-seq

- Melanoma cell line SKmel147
- Genes only from chromosome 5
- CRISPR ARID2KO vs ARID2 WT cell lines

### Deseq results table

File name	Cell line
Deseq_Results_ARID2Ko_v_ARID2wt_chr5.csv	SKmel147 (Melanoma)

# Script preparation

```
#define project directory
projectdir="/home/training/project_dir"

#practical data dir
practical_dir="/home/training/bulk_epigenetics_practical"

#container directory
container_dir=$practical_dir/"container/"

#blacklist_dir=
bl_regions=$practical_dir/"atacseq_data/supporting_files/hg38-blacklist.v2.bed"

#aligned_data
bowtie2_dir=$practical_dir/"atacseq_data/bowtie2"

# create directories for each packge in the data directory before running script
samtools_dir=$projectdir/"samtools"
picard_dir=$projectdir/"picard"
macs_dir=$projectdir/"macs"
deeptools_dir=$projectdir/"deeptools"
subread_dir=$projectdir/"subread"
homer_dir=$projectdir/"homer"
```

# Alignment: Bowtie2

- Alignment to GRCH38 with bowtie2

```
bowtie2 -p 12 \  
-x /sc/arion/projects/BiNGS/bings_omics/engine/annotation/homo_sapiens/grch38_gencode_36/bowtie2/2.2.8/index/index \  
-1 ${trimmgalore_dir}/${samplename}_ATAC_chr5_R1_val_1.fq.gz \  
-2 ${trimmgalore_dir}/${samplename}_ATAC_chr5_R2_val_2.fq.gz \  
-S ${bowtie2_dir}/${samplename}_noA.sam \  
-X 2000 \  
2> ${bowtie2_dir}/${samplename}.log |
```



# Alignment: Bowtie2

## ARID2 WT

```
5407385 reads; of these:
  5407385 (100.00%) were paired; of these:
    79578 (1.47%) aligned concordantly 0 times
    3250212 (60.11%) aligned concordantly exactly 1 time
    2077595 (38.42%) aligned concordantly >1 times
    ----
    79578 pairs aligned concordantly 0 times; of these:
      12484 (15.69%) aligned discordantly 1 time
    ----
    67094 pairs aligned 0 times concordantly or discordantly; of these:
      134188 mates make up the pairs; of these:
        82711 (61.64%) aligned 0 times
        28632 (21.34%) aligned exactly 1 time
        22845 (17.02%) aligned >1 times
99.24% overall alignment rate
```

## ARID2 KO

```
5722715 reads; of these:
  5722715 (100.00%) were paired; of these:
    85307 (1.49%) aligned concordantly 0 times
    3211482 (56.12%) aligned concordantly exactly 1 time
    2425926 (42.39%) aligned concordantly >1 times
    ----
    85307 pairs aligned concordantly 0 times; of these:
      14640 (17.16%) aligned discordantly 1 time
    ----
    70667 pairs aligned 0 times concordantly or discordantly; of these:
      141334 mates make up the pairs; of these:
        88304 (62.48%) aligned 0 times
        29711 (21.02%) aligned exactly 1 time
        23319 (16.50%) aligned >1 times
99.23% overall alignment rate
```

# Filtering/Sorting/Indexing and removing Dups



- Filtering/Sorting/Indexing

```
mkdir -p $samtools_dir
start_time=`date +%s`

#remove chrM and keep only "chr" chromosomes
singularity exec $container_dir/"samtools v1.15.sif" \
samtools view -h ${bowtie2_dir}/${samplename}_noA.sam | awk '($3 != "chrM" && $3 != "chrUn")' > ${bowtie2_dir}/${samplename}_noA_noM_chr_q20.sam

#convert to bam
singularity exec $container_dir/"samtools v1.15.sif" \
samtools view -h -q 20 ${bowtie2_dir}/${samplename}_noA_noM_chr_q20.sam -o ${samtools_dir}/${samplename}_noA_noM_chr_q20.bam

#sort
singularity exec $container_dir/"samtools v1.15.sif" \
samtools sort -@ 8 ${samtools_dir}/${samplename}_noA_noM_chr_q20.bam -o ${samtools_dir}/${samplename}_noA_noM_chr_q20_sorted.bam

#index
singularity exec $container_dir/"samtools v1.15.sif" \
samtools index -@ 8 -b ${samtools_dir}/${samplename}_noA_noM_chr_q20_sorted.bam

end_time=`date +%s`
echo execution time was `expr $end_time - $start_time` s.
```

- Removing duplicates

```
82 mkdir -p $picard_dir
83 start_time=`date +%s`
84
85 singularity exec ${container_dir}/picard_v2.9.2.sif \
86 java -jar /usr/local/share/picard-2.9.2-1/picard.jar MarkDuplicates \
87     I=${samtools_dir}/${samplename}_noA_noM_chr_q20_sorted.bam \
88     O=${picard_dir}/${samplename}_noA_noM_chr_q20_sorted_nd.bam \
89     REMOVE_DUPLICATES=true \
90     VALIDATION_STRINGENCY=LENIENT \
91     M=${picard_dir}/${samplename}_noA_noM_chr_q20_sorted_nd.txt
92
93 end_time=`date +%s`
94 echo execution time was `expr $end_time - $start_time` s.
95
```

# Filtering/Sorting/Indexing and removing Dups



- Re-sort and re-index

```
start_time=`date +%s`

singularity exec $container_dir/"samtools_v1.15.sif" \
samtools sort -@ 8 ${picard_dir}/${samplename}_noA_noM_chr_q20_sorted_nd.bam -o ${samtools_dir}/${samplename}_final.bam

singularity exec $container_dir/"samtools_v1.15.sif" \
samtools index -@ 8 -b ${samtools_dir}/${samplename}_final.bam

end_time=`date +%s`
echo execution time was `expr $end_time - $start_time` s.
```

- Generate bigwigs

```
mkdir -p $deeptools_dir
start_time=`date +%s`

singularity exec $container_dir/"deeptools_v3.5.1.sif" \
bamCoverage --bam ${samtools_dir}/${samplename}_final.bam \
             --outFileName $deeptools_dir/${samplename}_final.bw \
             --outFileFormat bigwig \
             --binSize=10 \
             --normalizeUsing RPKM \
             --extendReads=200 \
             --numberOfProcessors 8

end_time=`date +%s`
echo execution time was `expr $end_time - $start_time` s.
```

# Merging Bams and calling peaks

- Merging Bam files

```
#####  
### Merging bam files#####  
#####  
  
echo "Merging bam files"  
  
singularity exec ${container_dir}/samtools_v1.15.sif \  
  samtools merge ${samtools_dir}/master_atac.bam ${samtools_dir}/*_final.bam  
  
singularity exec ${container_dir}/samtools_v1.15.sif \  
  samtools index -@ 12 -b ${samtools_dir}/master_atac.bam
```

- Calling peaks

```
singularity exec ${container_dir}/macs_v2.2.9.1.sif \  
  macs2 callpeak --nomodel \  
    -t ${samtools_dir}/master_atac.bam \  
    --outdir ${macs_dir} \  
    -n master_atac \  
    -f BAMPE \  
    -g hs \  
    --keep-dup all \  
    --slocal 1000 \  
  2> ${macs_dir}/master_atac_mac2.log  
  
end time=`date +%s`
```

- Master bam peaks = 8026 (after blacklisting)

# Peak filtering and coverage tracks



- Peak filtering/formatting

```
#####
# remove bl regions and generating bed file for gb
#####

echo "Peak filtering/formatting"

singularity exec ${container_dir}/bedtools_v2.31.1.sif \
  bedtools intersect -a ${macs_dir}/master_atac_peaks.narrowPeak -b $bl_regions -v > ${macs_dir}/master_atac_peaks_bl.narrowPeak

awk '{print $1,$2,$3,$4}' ${macs_dir}/master_atac_peaks_bl.narrowPeak > ${macs_dir}/temp.bed

echo "track name=\"SKmell47_ATAC_master_regions\" description=\"SKmell47_ATAC_master_regions\"" | cat - ${macs_dir}/temp.bed > ${macs_dir}/master_atac_peaks_bl.bed
rm ${macs_dir}/temp.bed
```

- Generation of coverage tracks

```
#####
### run deeptools bamcoverage to generate bigwig
#####


echo "generating bigwig: master bam"

start_time=`date +%s`
singularity exec ${container_dir}/deeptools_v3.5.1.sif \
  bamCoverage --bam ${samtools_dir}/master_atac.bam \
               --outFileName $deeptools_dir/master_atac.bw \
               --outFileFormat bigwig \
               --binSize=10 \
               --normalizeUsing RPKM \
               --extendReads=200 \
               --numberOfProcessors 12

end_time=`date +%s`
echo execution time was `expr $end_time - $start_time` s.
```

# Coverage Track Visualization: IGV

- Adding bigwigs to IGV



# Quantification of RIPs

- Quantification of reads in peaks

```
# #####
# #### run featurecounts to generate coun matrix
# #####

echo "generating count matrix"
start_time=`date +%s`

awk 'OFS="\t" {print $4, $1, $2, $3, "."}' ${macs_dir}/master_atac_peaks_bl.narrowPeak > ${macs_dir}/master_atac_peaks_bl.saf

mkdir -p $subread_dir

singularity exec ${container_dir}/subread v2.0.1.sif \
featureCounts -p -a ${macs_dir}/master_atac_peaks_bl.saf -F SAF -o ${subread_dir}/master_atac_peaks_bl_subread.txt ${samtools_dir}/*final.bam -T 12
```


KO    WT

Assigned	977981	1161417
Unassigned_Unmapped	0	0
Unassigned_Read_Type	0	0
Unassigned_Singleton	0	0
Unassigned_MappingQuality	0	0
Unassigned_Chimera	0	0
Unassigned_FragmentLength	0	0
Unassigned_Duplicate	0	0
Unassigned_MultiMapping	0	0
Unassigned_Secondary	0	0
Unassigned_NonSplit	0	0
Unassigned_NoFeatures	951625	1050525
Unassigned_Overlapping_Length	0	0
Unassigned_Ambiguity	307	471



# Heatmap Visualization: Deeptools

```
#####  
### run deeptools to generate heatmap (split by sample)  
#####  
echo "generating deeptool heatmap - individual sample"  
  
start_time=`date +%s`  
  
awk 'FNR==NR {a[$4]; next} FNR> 1 && $4 in a' ${macs_dir}/master_atac_peaks_bl.bed ${macs_dir}/master_atac_summits.bed > ${macs_dir}/master_atac_summits_bl.bed  
  
singularity exec ${container_dir}/deeptools_v3.5.1.sif \  
computeMatrix reference-point \  
-R ${macs_dir}/master_atac_summits_bl.bed \  
--skipZeros \  
-S ${deeptools_dir}/SKMel147_ARID2WT_final.bw ${deeptools_dir}/SKMel147_ARID2KO_final.bw \  
-o ${deeptools_dir}/master_atac_summits_bl_2_samples.gz \  
-b 5000 -a 5000 \  
--referencePoint center \  
-p 4 \  
--samplesLabel ARID2WT ARID2KO  
  
singularity exec ${container_dir}/deeptools_v3.5.1.sif \  
plotHeatmap -m ${deeptools_dir}/master_atac_summits_bl_2_samples.gz \  
--plotFileFormat pdf \  
-out ${deeptools_dir}/master_atac_summits_bl_2_samples.pdf \  
--outFileSortedRegions ${deeptools_dir}/master_atac_summits_bl_2_samples_sorted.bed \  
--dpi 720 \  
--missingDataColor White \  
--colorMap Reds \  
--regionsLabel peaks \  
--heatmapHeight 13
```



# Heatmap Visualization: Deeptools

```
#####  
### run deeptools to generate heatmap (merged)  
#####  
  
echo "generating deeptool heatmap - merged sample"  
  
start_time=`date +%s`  
  
awk 'FNR==NR {a[$4]; next} FNR> 1 && $4 in a' ${macs_dir}/master_atac_peaks_bl.bed ${macs_dir}/master_atac_summits.bed > ${macs_dir}/master_atac_summits_bl.bed  
  
singularity exec ${container_dir}/deeptools_v3.5.1.sif \  
computeMatrix reference-point \  
-R ${macs_dir}/master_atac_summits_bl.bed \  
--skipZeros \  
-S ${deeptools_dir}/master_atac.bw \  
-o ${deeptools_dir}/master_atac_summits_bl_merged.gz \  
-b 5000 -a 5000 \  
--referencePoint center \  
-p 4 \  
--samplesLabel master_atac  
  
singularity exec ${container_dir}/deeptools_v3.5.1.sif \  
plotHeatmap -m ${deeptools_dir}/master_atac_summits_bl_merged.gz \  
--plotFileFormat pdf \  
-out ${deeptools_dir}/master_atac_summits_bl_merged.pdf \  
--outFileSortedRegions ${deeptools_dir}/master_atac_summits_bl_merged_sorted.bed \  
--dpi 720 \  
--missingDataColor White \  
--colorMap Reds \  
--regionsLabel peaks \  
--heatmapHeight 13
```

# Normalization and differential peak analysis



- With AWK

```
#####  
### Normalize and compute log2(KO/WT) add differential status  
#####  
  
tail -n +3 ${subread_dir}/master_atac_peaks_bl_subread.txt > ${subread_dir}/master_atac_peaks_bl_subread_no_header.txt  
  
#calculate total reads per sample  
total_reads_wt=$(awk '{ sum += $8 } END { print sum }' ${subread_dir}/master_atac_peaks_bl_subread_no_header.txt)  
total_reads_ko=$(awk '{ sum += $7 } END { print sum }' ${subread_dir}/master_atac_peaks_bl_subread_no_header.txt)  
  
#replace 0s  
awk '{ if ($7 == 0) $7 = 1; print }' ${subread_dir}/master_atac_peaks_bl_subread_no_header.txt > ${subread_dir}/master_atac_peaks_bl_subread_no_header_nozero.txt  
  
#calculate rpkm for wt  
awk -v total_reads_wt=$total_reads_wt '{ printf "%s\t%s\t%.2f\t%.2f\n", $1, $6, ($8 / ($6 / 1000)) / (total_reads_wt / 1000000), $7 }' ${subread_dir}/master_atac_peaks_bl_subread_no_header_nozero.txt > ${subread_dir}/master_atac_peaks_bl_subread_norm.txt  
awk -v total_reads_ko=$total_reads_ko '{ printf "%s\t%s\t%.2f\t%.2f\n", $1, $2, $3, ($4 / ($2 / 1000)) / (total_reads_ko / 1000000) }' ${subread_dir}/master_atac_peaks_bl_subread_norm.txt > ${subread_dir}/master_atac_peaks_bl_subread_norm_final.txt  
  
#calculate log2fc  
awk '{ ($5 = log($4/$3)/log(2)) print }' ${subread_dir}/master_atac_peaks_bl_subread_norm_final.txt > ${subread_dir}/master_atac_peaks_bl_subread_no_header_log2diff.txt  
  
# save up in KO peaks  
awk '$5 >= 1' ${subread_dir}/master_atac_peaks_bl_subread_no_header_log2diff.txt > ${subread_dir}/master_atac_peaks_bl_up_KO.txt  
awk 'FNR==NR {a[$1]; next} FNR> 1 && $4 in a' ${subread_dir}/master_atac_peaks_bl_up_KO.txt ${macs_dir}/master_atac_summits.bed > ${subread_dir}/master_atac_peaks_bl_up_KO_summit.bed  
  
# save down in KO peaks  
awk '$5 <= -1' ${subread_dir}/master_atac_peaks_bl_subread_no_header_log2diff.txt > ${subread_dir}/master_atac_peaks_bl_down_KO.txt  
awk 'FNR==NR {a[$1]; next} FNR> 1 && $4 in a' ${subread_dir}/master_atac_peaks_bl_down_KO.txt ${macs_dir}/master_atac_summits.bed > ${subread_dir}/master_atac_peaks_bl_down_KO_summit.bed  
  
# save static peaks  
awk '$5 < 1 && $5 > -1' ${subread_dir}/master_atac_peaks_bl_subread_no_header_log2diff.txt > ${subread_dir}/master_atac_peaks_bl_static.txt  
awk 'FNR==NR {a[$1]; next} FNR> 1 && $4 in a' ${subread_dir}/master_atac_peaks_bl_static.txt ${macs_dir}/master_atac_summits.bed > ${subread_dir}/master_atac_peaks_bl_static_summit.bed
```

# Visualization: Differential peak analysis

```
#####
### plot heatmap on differential and static regions
#####

echo "generating deeptool heatmap - differential and static peaks"

start_time=`date +%s`

for file in ${subread_dir}/*_summit.bed;
do
# file="$(basename $file)"
# samplename=$(echo "$file" | sed 's/.*\(_bl_\)/\1/g')
samplename=$(echo "$file" | sed 's/.*_bl_//g')

singularity exec ${container_dir}/deeptools_v3.5.1.sif \
computeMatrix reference-point \
-R $file \
--skipZeros \
-S ${deeptools_dir}/SKMell147_ARID2WT_final.bw ${deeptools_dir}/SKMell147_ARID2KO_final.bw \
-o ${deeptools_dir}/${samplename}.gz \
-b 5000 -a 5000 \
--referencePoint center \
-p 4 \
--samplesLabel ARID2WT ARID2KO

singularity exec ${container_dir}/deeptools_v3.5.1.sif \
plotHeatmap -m ${deeptools_dir}/${samplename}.gz \
--plotFileFormat pdf \
-out ${deeptools_dir}/${samplename}.pdf \
--dpi 720 \
--missingDataColor White \
--colorMap Reds \
--regionsLabel ${samplename} \
--heatmapHeight 13


done

end_time=`date +%s`
echo execution time was `expr $end_time - $start_time` s.
```




# Visualization: Differential peak analysis


N = 762



N = 496



N = 6765







# Motif analysis: Differential peak analysis





```
#####  
### homer motif analysis - split by differential and static peaks  
#####  
  
echo "motif enrichment analysis using homer"  
  
preparedDir=${practical_dir}/"atacseq_data/supporting_files/prepared_dir"  
  
start_time=`date +%s`  
  
for file in ${subread_dir}/*_summit.bed;  
do  
    samplename=$(echo "$file" | sed 's/.*bl.//g' | sed 's/_summit.bed//')  
    output_dir=${homer_dir}/${samplename}  
    mkdir -p $output_dir  
  
    bed_path=$file  
  
    singularity exec ${container_dir}/homer_v4.11_hg38.sif \\  
    findMotifsGenome.pl $bed_path hg38 $output_dir -preparedDir $preparedDir -size 200 -p 6  
done
```

# Motif analysis: Differential peak analysis

- Homer de-novo UP in KO

Rank	Motif	P-value	log P-pvalue	% of Targets	% of Background	STD(Bg STD)	Best Match/Details	Motif File
1		1e-40	-9.412e+01	29.32%	4.39%	47.9bp (65.9bp)	BATF(bZIP)/Th17-BATF-ChIP-Seq(GSE39756)/Homer(0.975) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
2		1e-26	-6.212e+01	14.66%	1.37%	42.5bp (58.7bp)	BORIS(Zf)/K562-CTCF-ChIP-Seq(GSE32465)/Homer(0.881) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
3		1e-12	-2.870e+01	16.54%	4.65%	56.5bp (62.8bp)	PB0086.1_Tcfap2b_1/Jaspar(0.722) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
4		1e-12	-2.834e+01	2.26%	0.01%	52.9bp (28.7bp)	PB0051.1_Osr2_1/Jaspar(0.646) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>

- Homer de-novo Down in KO

Rank	Motif	P-value	log P-pvalue	% of Targets	% of Background	STD(Bg STD)	Best Match/Details	Motif File
1		1e-211	-4.880e+02	36.88%	5.41%	43.7bp (65.8bp)	Fra1(bZIP)/BT549-Fra1-ChIP-Seq(GSE46166)/Homer(0.993) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
2		1e-29	-6.809e+01	7.19%	1.41%	51.8bp (59.4bp)	RUNX-AML(Runt)/CD4+ -PolII-ChIP-Seq(Barski_et_al.)/Homer(0.943) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
3		1e-25	-5.804e+01	38.10%	23.70%	52.2bp (62.8bp)	EWS:ERG-fusion(ETS)/CADO_ES1-EWS:ERG-ChIP-Seq(SRA014231)/Homer(0.945) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>
4		1e-18	-4.249e+01	3.45%	0.51%	51.0bp (61.9bp)	IRF1(IRF)/PBMC-IRF1-ChIP-Seq(GSE43036)/Homer(0.950) <a href="#">More Information</a>   <a href="#">Similar Motifs Found</a>	<a href="#">motif file (matrix)</a>



# Peak annotation

```
#####  
### associate peaks with nearest gene promoter  
#####  
  
##Peak annotation  
echo "annotate peaks with nearest gene and genomic localization"  
  
for file in ${subread_dir}/*_summit.bed;  
  
do  
    samplename=$(echo "$file" | sed 's/.*bl.//g' | sed 's/_summit.bed//')  
  
    singularity exec ${container_dir}/homer_v4.11_hg38.sif \  
    annotatePeaks.pl ${file} hg38 -size 200 > ${subread_dir}/${samplename}_annotated.txt  
  
done
```

# Pathway analysis

<https://maayanlab.cloud/Enrichr/>



## Enrichr

[Login](#) | [Register](#)

65,207,211 sets analyzed

494,081 terms

225 libraries

Analyze

[What's new?](#)

[Libraries](#)

[Gene search](#)

[Term search](#)

[About](#)

[Help](#)

## Input data

Expand a gene, a term, or a variant into a gene set:

e.g. STAT3, breast cancer, or rs28897756



Try an example [STAT3](#) [breast cancer](#) [rs28897756](#)

Include the top 100 most relevant genes

Paste a set of Entrez gene symbols on each row in the textbox below. You can try a gene set **example**. Also, you can now try adding a **background**.

Paste a set of valid Entrez gene symbols (e.g. STAT3) on each row in the text-box

0 gene(s) entered

In order to enable others to search your set please enter a brief description of it.

☐ Contribute your set so it can be searched by others

Submit

## Down in KO

GO Biological Process 2023

Bar Graph

Table

Clustergram

Appyter

Click the bars to sort. Now sorted by p-value ranking.

SVG PNG JPG

Taurine Metabolic Process (GO:0019530)

Mesonephros Development (GO:0001823)

Exonucleolytic Catabolism Of Deadenylated mRNA (GO:0043928)

Adenylate Cyclase-Activating Adrenergic Receptor Signaling Pathway (GO:0071880)

Positive Regulation Of Receptor Signaling Pathway Via STAT (GO:1904894)

Regulation Of Neurotransmitter Secretion (GO:0046928)

Calcium-Dependent Cell-Cell Adhesion Via Plasma Membrane Cell Adhesion Molecules (GO:0016339)

B Cell Activation Involved In Immune Response (GO:0002312)

Nuclear-Transcribed mRNA Catabolic Process, Exonucleolytic (GO:0000291)

Adrenergic Receptor Signaling Pathway (GO:0071875)

## Up in KO

GO Biological Process 2023

Bar Graph

Table

Clustergram

Appyter

Click the bars to sort. Now sorted by p-value ranking.

SVG PNG JPG

Calcium-Dependent Cell-Cell Adhesion Via Plasma Membrane Cell Adhesion Molecules (GO:0016339)

Embryo Development Ending In Birth Or Egg Hatching (GO:0009792)

Cell-Cell Adhesion Mediated By Cadherin (GO:0044331)

Adherens Junction Organization (GO:0034332)

Protein Localization To Site Of Double-Strand Break (GO:1990166)

TORC2 Signaling (GO:0038203)

Insulin Processing (GO:0030070)

Negative Regulation Of Cell-Substrate Adhesion (GO:0010812)

Macropinocytosis (GO:0044351)

Positive Regulation Of Mitophagy (GO:1901526)

# Heatmap with ChIPs - Kmeans

